



FRONT VIEW



SIDE VIEW



BOTTOM VIEW

FIG. 20.—Mechanical drawing of Cabot posterior wire leg splint. To be used with or without side splints. Uses: Injuries to the soft parts of the lower limb requiring fixation in transport; slight injuries to the knee or ankle requiring fixation in transport; fractures of the fibula; wounds of the ankle joint; injuries to the foot. The rods of the splint should be thoroughly padded, and they may be bent to allow flexion at the knee if desired. Side splint of wood or wire ladder may be used in connection with the splint if desired.