



FIG. 154.—EXTRACTION OF A FOREIGN BODY UNDER FLUOROSCOPIC CONTROL. THE OPEN SCREEN METHOD IN DARKENED ROOM. THE RED LIGHT ABOVE THE TABLE SERVES TO ILLUMINATE THE SURGICAL FIELD IN THE PRELIMINARY STEPS BEFORE THE DIRECT SEARCH FOR THE FOREIGN BODY BEGINS. BY MEANS OF THE FOOT SWITCH THE RED LIGHT IS TURNED OUT AND CURRENT IS TURNED INTO THE X-RAY TUBE BENEATH THE TABLE, AS ILLUSTRATED IN FIGURE 155